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CE
1639

Eakin Surgical Ltd, Greypoint, Cardiff Business Park,
Cardiff CF14 5WF, UK
+44 (0) 2920 767 800 info@eakinsurgical.com

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Products covered:

Kerrison Punches - S125 Kerrison Punch 1mm x 200mm 45°; S126 Kerrison Punch 2mm x 200mm 45°; S127 Kerrison Punch 3mm x 200mm 45°; S137 Kerrison Punch 4mm x 200mm 45°; S138 Kerrison Punch 5mm x 200mm 45°; S161 Kerrison Punch 1mm x 200mm 45° for Delicate Use; S261 Kerrison Punch 2mm x 200mm 45° for Delicate Use; S262 Kerrison Punch 3mm x 200mm 45° for Delicate Use; S367 Kerrison Punch 2mm x 170mm 90°.

Biopsy Punches - S263 Tischler-Morgan Biopsy Punch 200mm Straight Tip; S320 Tischler-Morgan Biopsy Punch 200mm Tip Down.

Symbol key:

| | | | |
|---|---|---|---|
|  | Manufacturer |  | Caution |
|  | Catalogue number |  | Medical device |
|  | Use-by-date |  | Unique device identifier |
|  | Date of manufacture |  | Do not resterilize |
|  | Authorised representative in the european community/ european union |  | CAUTION: US Federal law restricts this device to sale by or on the order of the physician |
|  | Batch code | QTY | Quantity of content |
|  | Single sterile barrier system with protective packaging |  | Indicates recovery/ recyclable package or package material |
|  | Do not reuse |  | Keep dry |
|  | Do not use if packaging is damaged and consult instructions for use |  | Consult instructions for use or consult electronic instructions for use |
|  | Sterilised using ethylene oxide |  | Single sterile barrier system |
|  | Latex Free |  | CE certification mark |

Intended Use:

Punch - A double action instrument used for cutting or extracting tough tissue, particularly bone. For use in general populations.

Contraindications:

- Not for use within the central nervous system or central circulatory system.

Intended user qualification and environment of use:

For use by trained medical professionals in a surgical setting with appropriate equipment.

⚠ Warnings:

- Inspect all sterile packaging prior to use and dispose of device if sterile barrier is broken or damaged.
- Do not use expired product.
- Do not use the device if damaged.
- Device is for single use only and cannot be reprocessed due to risks associated with infection and cross contamination. Cutting edges are unlikely to perform as effectively on re-use.
- Use with care in conjunction with heat sources or lasers to avoid injury to surgical staff or patient.
- Not for use in conjunction with magnetic therapies (MRI).

Use instructions:

- Inspect all sterile packaging prior to use and dispose of device if sterile barrier is broken or damaged.
- Inspect label for shelf-life information.
- Remove the device in accordance with sterile practise.
- No assembly of the device is required.
- Inspect device prior to use, damaged devices should be reported to the manufacturer and not used. Unused damaged items should be returned to the manufacturer.
- Following completion of the procedure discard the device in line with hospital guidelines.

Storage and Disposal:

- Store in a dry environment within shelf life as marked. Do not open until required.
- Dispose of the device as clinical waste and treat as sharps in accordance with hospital guidelines.

Manufacturer Details:

Contact the manufacturer for technical support. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority

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Eakin B.V. Medelse Poort, Biezenwei 23, 4004 MB,