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CE
1639

Eakin Surgical Ltd, Greypoint, Cardiff Business Park,
Cardiff CF14 5WF, UK
+44 (0) 2920 767 800 info@eakinsurgical.com

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
Intended Use:

Single use kit for myringotomy procedures.

Symbol key:

	Manufacturer		Caution
	Catalogue number		Medical device
	Use-by-date		Unique device identifier
	Date of manufacture		Do not re sterilize
	Authorised representative in the european community/ european union		CAUTION: US Federal law restricts this device to sale by or on the order of the physician
	Batch code	QTY	Quantity of content
	Single sterile barrier system with protective packaging		Indicates recovery/ recyclable package or package material
	Do not reuse		Keep dry
	Do not use if packaging is damaged and consult instructions for use		Consult instructions for use or consult electronic instructions for use
	Sterilised using ethylene oxide		Single sterile barrier system
	Latex Free		CE certification mark

Contraindications:

- Not for use within the central nervous system or central circulatory system.

Intended user qualification and environment of use:

For use by trained medical professionals in a surgical setting with appropriate equipment.

⚠ Warnings:

- Inspect all sterile packaging prior to use and dispose of device if sterile barrier is broken or damaged.
- Do not use expired product.
- Do not use the device if damaged.
- Device is for single use only and cannot be reprocessed due to risks associated with infection and cross contamination. Cutting edges are unlikely to perform as effectively on re-use.
- Use with care in conjunction with heat sources or lasers to avoid injury to surgical staff or patient.
- Not for use in conjunction with magnetic therapies (MRI).

Use instructions:

- Inspect all sterile packaging prior to use and dispose of device if sterile barrier is broken or damaged.
- Inspect label for shelf-life information.
- Remove the device in accordance with sterile practise.
- No assembly of the device is required.
- Inspect device prior to use, damaged devices should be reported to the manufacturer and not used. Unused damaged items should be returned to the manufacturer.
- Following completion of the procedure discard the device in line with hospital guidelines.

Storage and Disposal:

- Store in a dry environment within shelf life as marked. Do not open until required.
- Dispose of the device as clinical waste and treat as sharps in accordance with hospital guidelines.

Manufacturer Details:

Contact the manufacturer for technical support. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority

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Eakin B.V. Medelse Poort, Biezenwei 23, 4004 MB,